



Metafix

PRODUCT INFORMATION BULLETIN

Hygea 2 HTM 01-05 Compliant Ultrasonic Cleaner



FEATURES AND BENEFITS

- ◆ 6–10 minute complete cycle
- ◆ Up to 45 instruments per cycle in 8L tank
- ◆ Uses 90%+ less energy than washer disinfectors
- ◆ Anti-microbial surface technology
- ◆ Menu Driven LCD control with admin/operator logins
- ◆ Auto locking door
- ◆ Micro SD Card Port for data logging and integral printer
- ◆ Approved on the NHS Supply Chain framework agreement
- ◆ Approved on the Scottish Procurement contract
- ◆ 3 year warranty
- ◆ Nationwide Metafix engineer support and service

Metafix (UK) Limited

1 Dean Close, Raunds, Northamptonshire NN9 6BD

E-mail: sales@metafix.co.uk

Tel: 01933 461907

www.metafix.co.uk

Hygea 2 HTM 01–05 Compliant Ultrasonic Cleaner



CLEANING CYCLES

- ◆ 6 – 10 minutes programmable
- ◆ 20°C to 40°C programmable

INCLUDED WITH ULTRASONIC CLEANER

- ◆ 3 tier basket tray system
- ◆ 3 year warranty

ACCESSORIES

- ◆ Hygea Valet instrument cassettes
- ◆ Sonozyme 3 triple enzymatic instrument cleaner sachets
- ◆ Sonozyme 3 triple enzymatic instrument cleaner 1L
- ◆ Aluminium Foil
- ◆ Hygea ultrasonic activity meter
- ◆ Load check indicators
- ◆ Residual Protein Test
- ◆ Printer paper rolls



SPECIFICATIONS

- ◆ Unit size (OD) 460(W) x 430(D) x 375(H) mm (lid open height 655mm)
- ◆ Tank size (ID) 285(W) x 205(D) x 160(H) mm
- ◆ Chamber capacity 8 litres/up to 45 standard dental instruments
- ◆ Instrument tray size Up to 3 trays each: 240(W) x 170(D) x 15(H) mm
- ◆ Electrical 220–240V, 200W operating
- ◆ Cleaning temperature 20°C to 40°C programmable
- ◆ Draining Easy draining tank – tap and hose
- ◆ Manufactured by Ultrawave UK to ISO9001/ISO13485
- ◆ Approvals CE HTM01–05 and HTM2030 compliant

Metafix (UK) Limited

1 Dean Close, Raunds, Northamptonshire NN9 6BD
E-mail: sales@metafix.co.uk

Tel: 01933 461907
www.metafix.co.uk